

## VDAC3 antibody

### Product Information

Catalog No.:	FNab09389
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background


Forms a channel through the mitochondrial outer membrane that allows diffusion of small hydrophilic molecules.

### Immunogen information


Immunogen:	voltage-dependent anion channel 3
Synonyms:	HD VDAC3, hVDAC3, VDAC 3, VDAC3
Observed MW:	23-27 kDa,45 kDa
Uniprot ID :	Q9Y277

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis using FNab09389(VDAC3 antibody) at dilution of 1:50



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab09389(VDAC3 antibody) at dilution of 1:300