

anti-Vertnin antibody

Product Information

Catalog No.:	FNab09395
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Vertnin, also named as VRTN and C14orf115, belongs to the vertnin family. It is an essential factor for development of the embryo in a wide range of organisms.

Immunogen information

Immunogen:	chromosome 14 open reading frame 115
Synonyms:	C14orf115, FLJ10811, Vertnin, VRTN
Observed MW:	78 kDa
UniprotID :	Q9H8Y1

Application


Reactivity:	Human
Tested Application:	ELISA, IF, IHC, WB
Recommended dilution:	WB: 1:200-1:1000; IHC: 1:20-1:200; IF: 1:10-1:100
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human placenta tissue slide using FNab09395(C14orf115 Antibody) at dilution of 1:50



SGC-7901 cells were subjected to SDS PAGE followed by western blot with FNab09395(C14orf115 Antibody) at dilution of 1:300