

## anti- VPS45 antibody

### Product Information

Catalog No.:	FNab09443
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

May play a role in vesicle-mediated protein trafficking from the Golgi stack through the trans-Golgi network.

### Immunogen information

Immunogen:	vacuolar protein sorting 45 homolog(S. cerevisiae)
Synonyms:	H VPS45, H1, H1VPS45, hVps45, VPS45, VPS45A, VPS45B, VPS54A, VSP45, VSP45A
Observed MW:	65 kDa
Uniprot ID :	Q9NRW7

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, IF
Recommended dilution:	IHC: 1:20-1:200; IF: 1:20-1:200
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



Immunofluorescent analysis of HepG2 cells, using FNab09443(VPS45 antibody) at 1:50 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent DNA dye).



Immunohistochemical of paraffin-embedded human testis using FNab09443(VPS45 antibody) at dilution of 1:50