

anti- VSTM2L antibody

Product Information

Catalog No.:	FNab09457
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	V-set and transmembrane domain containing 2 like
Synonyms:	C20orf102, dJ1118M15.2, VSTM2L
Observed MW:	22 kDa
Uniprot ID :	Q96N03

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:5000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Recombinant proteins were subjected to SDS PAGE followed by western blot with FNab09457(VSTM2L Antibody) at dilution of 1:1000