
anti- WBP1 antibody

Product Information

Catalog No.:	FNab09472
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	WW domain binding protein 1
Synonyms:	WBP 1, WBP1, WW domain binding protein 1
Observed MW:	32 kDa
UniprotID :	Q96G27

Application


Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:WB:	1:200-1:1000,IHC: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human pancreas cancer using FNab09472(WBP1 antibody) at dilution of 1:50



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab09472(WBP1 antibody) at dilution of 1:400