

anti- WDR6 antibody

Product Information

Catalog No.:	FNab09499
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Enhances the STK11/LKB1-induced cell growth suppression activity. Negative regulator of amino acid starvation-induced autophagy.

Immunogen information

Immunogen:	WD repeat domain 6
Synonyms:	WDR6, FLJ56107, FLJ52552, WD repeat domain 6, WD repeat containing protein 6
Observed MW:	122 kDa
UniprotID :	Q9NNW5

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:200-1:2000; IF: 1:10-1:100
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



HeLa cells were subjected to SDS PAGE followed by western blot with FNab09499(WDR6 antibody) at dilution of 1:500