

anti- XIRP2 antibody

Product Information

| | |
|---------------|---|
| Catalog No.: | FNab09542 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| Clone ID: | None |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |

Background

Protects actin filaments from depolymerization.

Immunogen information

| | |
|--------------|--|
| Immunogen: | xin actin-binding repeat containing 2 |
| Synonyms: | XIRP2, Beta xin, FLJ39111, Xeplin, DKFZp779G1747 |
| Observed MW: | 280 kDa |
| UniprotID : | A4UGR9 |

Application


| | |
|-----------------------|-----------------------------------|
| Reactivity: | Human, Mouse |
| Tested Application: | ELISA, WB, IHC |
| Recommended dilution: | WB: 1:200-1:2000; IHC: 1:20-1:200 |
| Image: | |

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human skeletal muscle using FNab09542(XIRP2 antibody) at dilution of 1:50



COLO 320 cells were subjected to SDS PAGE followed by western blot with FNab09542(XIRP2 antibody) at dilution of 1:500