

anti- XPO6 antibody

Product Information

Catalog No.:	FNab09550
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Mediates the nuclear export of actin and profilin-actin complexes in somatic cells.

Immunogen information

Immunogen:	exportin 6
Synonyms:	EXP6, exportin 6, FLJ22519, KIAA0370, Ran binding protein 20, RANBP20, XPO6
Observed MW:	129 kDa
UniprotID :	Q96QU8

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, IP
Recommended dilution:	IHC: 1:10-1:100; IP: 1:200-1:1000
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



IP Result of anti-XPO6 (IP:FNab09550, 4ug;
Detection:FNab09550 1:300) with HeLa cells
lysate 4000ug.



Immunohistochemical of paraffin-embedded
human ovary using FNab09550(XPO6 antibody) at
dilution of 1:25