

anti- ZAP70 antibody

Product Information

Catalog No.:	FNab09584
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	7A7
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	zeta-chain(TCR) associated protein kinase 70kDa
Synonyms:	70 kDa zeta associated protein, FLJ17670, FLJ17679, SRK, STD, Syk related tyrosine kinase, Tyrosine protein kinase ZAP 70, TZK, ZAP 70, ZAP70
Observed MW:	66 kDa
Uniprot ID :	P43403

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:1000-1:5000; IHC: 1:200-1:1000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human lymphoma using FNab09584(ZAP70 antibody) at dilution of 1:200



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab09584(ZAP70 antibody) at dilution of 1:3000