

anti- ZBTB20 antibody

Product Information

Catalog No.:	FNab09589
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

May be a transcription factor that may be involved in hematopoiesis, oncogenesis, and immune responses.

Immunogen information

Immunogen:	zinc finger and BTB domain containing 20
Synonyms:	DKFZp566F123, DPZF, HOF, ODA 8S, ZBTB20, Zinc finger protein 288, ZNF288
Observed MW:	73 kDa
Uniprot ID :	Q9HC78

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:5000; IHC: 1:20-1:200
Image:	

1


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human kidney tissue slide using FNab09589(ZBTB20 Antibody) at dilution of 1:50



L02 cells were subjected to SDS PAGE followed by western blot with FNab09589(ZBTB20 antibody) at dilution of 1:1000