

anti- ZBTB46 antibody

Product Information

Catalog No.:	FNab09596
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Function as a transcriptional repressor for PRDM1.

Immunogen information

Immunogen:	zinc finger and BTB domain containing 46
Synonyms:	BTBD4, dJ583P15.7, dJ583P15.8, FLJ13502, Homo sapiens(Human), RINZF, ZBTB46, ZNF340
Observed MW:	65 kDa
UniprotID :	Q86UZ6

Application


Reactivity:	Human, Rat, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:1000
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



RAW 264.7 cells were subjected to SDS PAGE followed by western blot with FNab09596(ZBTB46 Antibody) at dilution of 1:600