

## anti- ZMPSTE24 antibody

### Product Information

Catalog No.:	FNab09649
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Proteolytically removes the C-terminal three residues of farnesylated proteins. Acts on lamin A/C.

### Immunogen information

Immunogen:	zinc metallopeptidase(STE24 homolog, <i>S. cerevisiae</i> )
Synonyms:	CAAX prenyl protease 1 homolog, FACE 1, FACE1, FLJ14968, STE24, Ste24p, ZMPSTE24
Observed MW:	50 kDa
UniprotID :	O75844

### Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:200-1:2000; IP: 1:200-1:1000; IF: 1:10-1:100
Image:	

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



PC-3 cells were subjected to SDS PAGE followed by western blot with FNab09649(ZMPSTE24 antibody) at dilution of 1:400

**Wuhan Fine Biotech Co., Ltd.**

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)