

anti- ZNF311 antibody

Product Information

Catalog No.:	FNab09683
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

May be involved in transcriptional regulation.

Immunogen information

Immunogen:	zinc finger protein 311
Synonyms:	Zf31, ZFP31, zinc finger protein 311, Zinc finger protein zfp 31, ZNF311
Observed MW:	76 kDa, observed 70 kDa and 60 kDa
UniprotID :	Q5JNZ3

Application


Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:1000-1:10000; IHC: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human testis slide using FNab09683(ZNF311 Antibody) at dilution of 1:50



HL-60 cells were subjected to SDS PAGE followed by western blot with FNab09683(ZNF311 antibody) at dilution of 1:4000