

anti- ZNF446 antibody

Product Information

Catalog No.:	FNab09700
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

May be involved in transcriptional regulation.

Immunogen information

Immunogen:	zinc finger protein 446
Synonyms:	FLJ20626, zinc finger protein 446, ZKSCAN20, ZNF446, ZSCAN30
Observed MW:	49 kDa
UniprotID :	Q9NWS9

Application


Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:200-1:1000; IF: 1:50-1:500
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with FNab09700(ZNF446 Antibody) at dilution of 1:300