

## anti- ZNF606 antibody

### Product Information

Catalog No.:	FNab09719
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

May act as a transcriptional repressor.

### Immunogen information

Immunogen:	zinc finger protein 606
Synonyms:	FLJ14260, KIAA1852, Zinc finger protein 328, zinc finger protein 606, ZNF328, ZNF606
Observed MW:	92 kDa
UniprotID :	Q8WXB4

### Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, IP, WB
Recommended dilution:	WB: 1:200-1:2000; IP: 1:200-1:1000
Image:	

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



IP Result of anti-ZNF606 (IP:FNab09719, 4ug;  
Detection:FNab09719 1:300) with mouse liver  
tissue lysate 4000ug.



mouse liver tissue were subjected to SDS PAGE  
followed by western blot with  
FNab09719(ZNF606 Antibody) at dilution of  
1:600