

anti- ZNF740 antibody

Product Information

Catalog No.:	FNab09735
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

May be involved in transcriptional regulation.

Immunogen information

Immunogen:	zinc finger protein 740
Synonyms:	TB7, Zfp740, zinc finger protein 740, ZNF740
Observed MW:	22 kDa
UniprotID :	Q8NDX6

Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, IP, WB
Recommended dilution:	WB: 1:200-1:2000; IP: 1:200-1:2000
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



IP Result of anti-ZNF740 (IP:FNab09735, 3ug;
Detection:FNab09735 1:600) with mouse brain
tissue lysate 4000ug.



mouse brain tissue were subjected to SDS PAGE
followed by western blot with
FNab09735(ZNF740 Antibody) at dilution of
1:600