

## IL2 antibody

### Product Information

Catalog No.:	FNab09820
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background


The protein encoded by this gene is a secreted cytokine that is important for the proliferation of T and B lymphocytes. The receptor of this cytokine is a heterotrimeric protein complex whose gamma chain is also shared by interleukin 4 (IL4) and interleukin 7 (IL7). The expression of this gene in mature thymocytes is monoallelic, which represents an unusual regulatory mode for controlling the precise expression of a single gene. The targeted disruption of a similar gene in mice leads to ulcerative colitis-like disease, which suggests an essential role of this gene in the immune response to antigenic stimuli.

### Immunogen information


Immunogen:	interleukin 2
Synonyms:	Aldesleukin, IL 2, IL2, interleukin 2, lymphokine, T cell growth factor, TCGF
Observed MW:	22 kDa
Uniprot ID :	P60568

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:100
Image:	



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using FNab09820( IL2 antibody) at dilution of 1:100



MCF-7 cells were subjected to SDS PAGE followed by western blot with FNab09820 at dilution of 1:1000