

anti- MEKK3 antibody

Product Information

Catalog No.:	FNab09830
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Component of a protein kinase signal transduction cascade. Mediates activation of the NF-kappa-B, AP1 and DDIT3 transcriptional regulators.

Immunogen information

Immunogen:	MEKK3
Synonyms:	MAP3K3, MAPK/ERK kinase kinase 3, MAPKKK3, MEK kinase 3, MEKK 3, MEKK3
Observed MW:	74 kDa
Uniprot ID :	Q99759

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:1000; IHC: 1:50-1:500
Image:	

1


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using FNab09830(MEKK3 antibody) at dilution of 1:200



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab09830(MEKK3 antibody) at dilution of 1:300