

anti- PDI antibody

Product Information

Catalog No.:	FNab09833	
Size:	100µg	
Form:	liquid	
Purification:	Protein A+G purification	
Purity:	\geq 95% as determined by SDS-PAGE	
Host:	Mouse	
Clonality:	monoclonal	
Clone ID:	3D9	
IsoType:	IgG2b	
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)	

Background

PDIA1(Protein disulfide-isomerase) is also named as ERBA2L, PDI, P4HB, PO4DB. It is a multifunctional protein that catalyzes the formation, breakage and rearrangement of disulfide bonds. In some cell types, it seems to be secreted or associated with the plasma membrane, where it undergoes constant shedding and replacement from intracellular sources. It can exsit as homodimer and monomers and homotetramers may also occur(PMID:12095988).

Immunogen information

Immunogen:	prolyl 4-hydroxylase, beta polypeptide
Synonyms:	DSI, ERBA2L, GIT, P4HB, P4Hbeta, p55, PDI, PDIA1, PHDB, PO4DB, PO4HB, PROHB, Protein disulfide isomerase
Observed MW:	58 kDa
Uniprot ID :	P07237

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889

www.fn-test.com



Recommended dilution: WB: 1:10000-1:100000; IHC: 1:100-1:400; IF: 1:50-1:500 Image:



Immunohistochemistry of paraffin-embedded human liver tissue slide using FNab09833(PDI Antibody) at dilution of 1:200

Immunofluorescent analysis of (-20°C Ethanol) fixed HeLa cells using FNab09833 (PDI antibody) at dilution of 1:100 and Alexa Fluor 488congugated Goat Anti-Mouse IgG(H+L)


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889

www.fn-test.com





HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab09833(PDI antibody) at dilution of 1:50000

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

www.fn-test.com