

anti- EPAC1 antibody

Product Information

Catalog No.:	FNab09888
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Rap guanine nucleotide exchange factor(GEF) 3
Synonyms:	RAPGEF3, CAMP-GEFI, EPAC, EPAC1, HSU79275, bcm910, rap guanine nucleotide exchange factor 3
Observed MW:	104 kDa
Uniprot ID :	O95398

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500 - 1:2000; IF: 1:50 - 1:200
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunofluorescence analysis of U2OS cells using FNab09888(EPAC1 antibody) at dilution of 1:100. Blue: DAPI for nuclear staining.



mouse heart were subjected to SDS PAGE followed by western blot with FNab09888(EPAC1 antibody) at dilution of 1:1000