

# anti- CTLA4 antibody

### **Product Information**

Catalog No.: FNab09937

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

## Immunogen information

Immunogen: cytotoxic T-lymphocyte protein 4

Synonyms: CTLA4, ALPS5, CD, CD152, CELIAC3, CTLA-4, GRD4, GSE,

IDDM12, cytotoxic T-lymphocyte protein 4

Observed MW: 37 kDa Uniprot ID: P16410

# **Application**

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, WB, IHC

Recommended dilution: WB: 1:500 - 1:2000; IHC: 1:20 - 1:200

Image:


1

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Immunohistochemistry of paraffin-embedded human tonsil tissue slide using FNab09937 (CTLA4 antibody) at dilution of 1:100

—100kd —70kd

-55kd

-40kd

—35kd

-25kd

-15kd

mouse lung were subjected to SDS PAGE followed by western blot with FNab09937 (CTLA4 antibody) at dilution of 1:1000