

anti- TREM1 antibody

Product Information

Catalog No.:	FNab09969
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Stimulates neutrophil and monocyte-mediated inflammatory responses. Triggers release of pro-inflammatory chemokines and cytokines, as well as increased surface expression of cell activation markers. Amplifier of inflammatory responses that are triggered by bacterial and fungal infections and is a crucial mediator of septic shock.

Immunogen information

Immunogen:	triggering receptor expressed on myeloid cells 1
Synonyms:	CD354,TREM-1,Triggering receptor expressed on monocytes 1
Observed MW:	30 kDa
Uniprot ID :	Q9NP99

Application

Reactivity:	Human, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com

Image:



HL-60 cells were subjected to SDS PAGE followed by western blot with FNab09969(TREM1 antibody) at dilution of 1:1000