

anti- TFAP2A antibody

Product Information

Catalog No.:	FNab09980
Size:	100µg
Form:	liquid
Purification:	protein A+G purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	6C5
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	transcription factor AP-2 alpha
Synonyms:	Activator protein 2, AP 2, AP 2 transcription factor, AP 2alpha, AP2, AP2 alpha, AP2TF, BOFS, FLJ51761, TFAP2, TFAP2A
Observed MW:	50 kDa
Uniprot ID :	P05549

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:2000 - 1:10000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



HeLa cells were subjected to SDS PAGE followed by western blot with FNab09980(TFAP2A antibody) at dilution of 1:5000