

# anti- HSPA5 antibody

## **Product Information**

Catalog No.:	FNab09987
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	$\geq$ 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### **Immunogen information**

Immunogen:	heat shock 70kDa protein 5
Synonyms:	BIP, GRP 78, GRP78, HSP A5, HSPA5, MIF2
Observed MW:	78 kDa
Uniprot ID :	P11021


# Application

Reactivity:	Human, Rat	
Tested Application:	ELISA, WB, IHC	
Recommended dilution: WB: 1:500-1:2000; IHC: 1:50-1:500		
Image:		


## Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)





Immunohistochemistry of paraffin-embedded human pancreas cancer using FNab09987(HSPA5 antibody) at dilution of 1:200



HeLa cells were subjected to SDS PAGE followed by western blot with FNab09987(HSPA5 antibody) at dilution of 1:1000

2

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)