

# anti- ABC3A antibody

## **Product Information**

Catalog No.:	FNab10020
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	$\geq$ 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### **Immunogen information**

Immunogen:	APOBEC-3A
Synonyms:	APOBEC3A, A3A, ARP3, PHRBN, bK150C2.1, DNA dC->dU-editing enzyme APOBEC-3A
Observed MW:	26 kDa
Uniprot ID :	P31941

## Application

Reactivity:	Human, Rat		
Tested Application:	ELISA, WB, IHC		
Recommended dilution: WB: 1:200-1:2000; IHC: 1:50-1:500			
Image:			


#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)


Tel :( 0086)027-87384275 Fax: (0086)027-87800889

www.fn-test.com





Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab10020(ABC3A Antibody) at dilution of 1:200



rat testis were subjected to SDS PAGE followed by western blot with FNab10020( ABC3A Antibody) at dilution of 1:1000

2

#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)