

# anti- CADH2 antibody

#### **Product Information**

Catalog No.: FNab10068

Size: 100µg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

# **Background**

This gene encodes a classical cadherin and member of the cadherin superfamily. Alternative splicing results in multiple transcript variants, at least one of which encodes a preproprotein is proteolytically processed to generate a calcium-dependent cell adhesion molecule and glycoprotein. This protein plays a role in the establishment of left-right asymmetry, development of the nervous system and the formation of cartilage and bone.

#### **Immunogen information**

Immunogen: cadherin 2, type 1, N-cadherin (neuronal)

Synonyms: Cadherin 2, CD325, CDH2, CDHN, CDw325, N cadherin, NCAD, N-

cadherin. Neural cadherin

Observed MW: 140 kDa Uniprot ID: P19022

## **Application**

Reactivity: Human, rat
Tested Application: ELISA, IHC

## Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)


Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>

1



Recommended dilution: IHC: 1:20-1:200

Image:



Immunohistochemistry of paraffin-embedded rat kidney using FNab10068( CADH2 Antibody) at dilution of 1:100