



anti- MYPR antibody

Product Information

Catalog No.:	FNab10080
Size:	100 μ g
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20 $^{\circ}$ C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Myelin proteolipid protein
Synonyms:	PLP1, GPM6C, HLD1, MMPL, PLP, PLP/DM20, PMD, SPG2; proteolipid protein 1
Calculated MW:	
Uniprot ID :	P60201

Application

Specificity:	Human, Rat
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



FineTest


Wuhan Fine Biotech Co.,Ltd.

C6-323 Biolake, No. 666 Gaoxin Ave.

East Lake High-tech Development District

Wuhan City, Hubei, China (Zip code 430072)

Tel: 0086-027-87384275 Fax:0086-027-87800889



Immunohistochemistry of paraffin-embedded rat brain tissue slide using FNab10080(MYPR Antibody) at dilution of 1:200