

anti- PTCH1 antibody

Product Information

Catalog No.: FNab10084

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

Acts as a receptor for sonic hedgehog(SHH), indian hedgehog(IHH) and desert hedgehog(DHH). Associates with the smoothened protein(SMO) to transduce the hedgehog's proteins signal. Seems to have a tumor suppressor function, as inactivation of this protein is probably a necessary, if not sufficient step for tumorigenesis.

Immunogen information

Immunogen: patched homolog 1

Synonyms: HPE7, NBCCS, Patched, patched homolog 1, Patched1, Protein patched

homolog 1, PTC, PTC1, PTCH, PTCH1, PTCH11

Observed MW: 161 kDa Uniprot ID: Q13635

Application

Reactivity: Human, Rat
Tested Application: ELISA, IHC

Recommended dilution: IHC: 1:20-1:200

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>

1



Image:



Immunohistochemistry of paraffin-embedded human colon cancer using FNab10084(PTCH1 antibody) at dilution of 1:50