

anti- FGR antibody

Product Information

Catalog No.:	FNab10137
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 24 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Tyrosine-protein kinase Fgr
Synonyms:	Gardner-Rasheed feline sarcoma viral (v-fgr) oncogene homolog, Proto- oncogene c-Fgr, p55-Fgr, p58-Fgr, FGR
Observed MW:	59 kDa
Uniprot ID :	P09769

Application


Reactivity:	Rat, Mouse; other species are not tested. Please decide the specificity by homology	
Tested Application:	ELISA, IHC	
Recommended dilution:IHC: 1:50-1:500		
Image:		

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

1





Immunohistochemistry of paraffin-embedded mouse spleen tissue slide using FNab10137(FGR Antibody) at dilution of 1:200

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

www.fn-test.com