

# anti- GRM2 antibody

#### **Product Information**

Catalog No.: FNab10138

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24

months (Avoid repeated freeze / thaw cycles.)

## Immunogen information

Immunogen: glutamate receptor, metabotropic 2

Synonyms: GLUR2, GPRC1B, GRM2, mGlu2, mGluR2

Observed MW: 96 kDa Uniprot ID: Q14416

## **Application**

Reactivity: Rat, Mouse; other species are not tested. Please decide the specificity by

homology

Tested Application: ELISA, IHC

Recommended dilution: IHC: 1:50-1:500

Image:


#### 1

## Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Immunohistochemistry of paraffin-embedded mouse brain tissue slide using FNab10138(GRM2 Antibody) at dilution of 1:200