

# anti- EWSR1 antibody

#### **Product Information**

Catalog No.: FNab10142

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

## **Background**

None

## **Immunogen information**

Immunogen: Ewing sarcoma breakpoint region 1

Synonyms: EWS, EWS oncogene, EWSR1, RNA binding protein EWS

Observed MW: 37kDa/68kDa/100kDa

Uniprot ID: Q01844

## **Application**

Reactivity: Human, Mouse Tested Application: ELISA, WB, IF

Recommended dilution: WB: 1:500-1:2000; IF: 1:50-1:200

Image:


#### 1

#### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Immunofluorescent analysis of U2OS cells using FNab10142 (EWSR1 antibody) at dilution of 1:100



mouse spleen were subjected to SDS PAGE followed by western blot with FNab10142(EWSR1 antibody) at dilution of 1:1000