

# anti- GSDMB antibody

#### **Product Information**

Catalog No.:	FNab10147
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	$\geq$ 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

The N-terminal moiety promotes pyroptosis. May be acting by homooligomerizing within the membrane and forming pores(PubMed:27281216). The physiological relevance of this observation is unknown(Probable).

### Immunogen information

Immunogen:	gasdermin B
Synonyms:	GSDMB, GSDML, PP4052, PRO2521, gasdermin-B
Observed MW:	47 kDa
Uniprot ID :	Q8TAX9

## Application

Reactivity:	Human, Mouse, Rat		
Tested Application:	WB		
Recommended dilution:WB: 1:200-1:2000			
Image:			


#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889

www.fn-test.com





Hela cells were subjected to SDS PAGE followed by western blot with FNab10147(GSDMB Antibody) at dilution of 1:1000

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

www.fn-test.com