

# anti- TLR8 antibody

# **Product Information**

Catalog No.:	FNab10155
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	$\geq$ 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

# Background

The protein encoded by this gene is a member of the Toll-like receptor (TLR) family which plays a fundamental role in pathogen recognition and activation of innate immunity. TLRs are highly conserved from Drosophila to humans and share structural and functional similarities. They recognize pathogen-associated molecular patterns (PAMPs) that are expressed on infectious agents, and mediate the production of cytokines necessary for the development of effective immunity. TLRs exhibit different patterns of expression. This gene is predominantly expressed in lung and peripheral blood leukocytes, and lies in close proximity to another family member, TLR7, on chromosome X.

## Immunogen information

Immunogen:	Toll-like receptor 8
Synonyms:	TLR8, CD288, TLR8
Observed MW:	120 kDa, 55 kDa
Uniprot ID :	Q9NR97


# Application

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)



Reactivity:Human, Mouse, RatTested Application:ELISA, WB, IHC, IFRecommended dilution:WB: 1:500 - 1:2000; IHC: 1:50 - 1:200; IF: 1:50 - 1:200Image:



Immunohistochemistry of paraffin-embedded human breast cancer tissue using FNab10155(TLR8 Antibody) at dilution of 1:100



A549 cells were subjected to SDS PAGE followed by western blot with FNab10155(TLR8 antibody) at dilution of 1:1000

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275

www.fn-test.com