

# anti-PIK3R2 antibody

#### **Product Information**

Catalog No.: FNab10160

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

## **Background**

Phosphatidylinositol 3-kinase (PI3K) plays an important role in the metabolic actions of insulin and is required for adipogenesis. The PI3K pathway has also been identified as an important player in cancer development and progression. The Class IA PI3K heterodimer, which is composed of a P110 catalytic subunit and a P85 regulatory subunit, is activated upon association of the P85 subunit with upstream adaptor proteins or receptor tyrosine kinases. Mutations in PIK3R1 are implicated in cases of breast cancer and associated to SHORT syndrome.

#### Immunogen information

Immunogen: Phosphatidylinositol 3-kinase regulatory subunit beta

Synonyms: MPPH, MPPH1, P85B, p85, p85-BETA, PIK3R2, PI3-kinase regulatory

subunit beta, PI3K regulatory subunit beta, PI3-kinase subunit p85-beta,

PtdIns-3-kinase regulatory subunit p85-beta

Observed MW: 81 kDa, 100 kDa

Uniprot ID: O00459

## **Application**

1

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)


Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>



Reactivity: Human, rat
Tested Application: ELISA, WB

Recommended dilution: WB: 1:500-1:2000

Image:



Hela cells were subjected to SDS PAGE followed by western blot with FNab10160(PIK3R2 antibody) at dilution of 1:1000