

anti- NOS3 antibody

Product Information

Catalog No.:	FNab10162
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	Nitric oxide synthase, endothelial
Synonyms:	Nanos homolog 3, nanos homolog 3 (Drosophila), NANOS1L, NANOS3, NOS 3, NOS3
Observed MW:	70 kDa, 133 kDa
Uniprot ID :	P29474

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:500
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded Rat kidney tissue slide using FNab10162(NOS3 antibody) at dilution of 1:200



239T cells were subjected to SDS PAGE followed by western blot with FNab10162(NOS3 antibody) at dilution of 1:1000