

anti- SNF5 antibody

Product Information

Catalog No.: FNab10186

Size: $100\mu g$ Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen: SWI/SNF-related matrix-associated actin-dependent regulator of

chromatin subfamily B member 1

Synonyms: BAF47, BRG1 associated factor 47, hSNF5, hSNFS, INI1, Integrase

interactor 1 protein, RDT, Sfh1p, SMARCB1, SNF5, SNF5 homolog,

SNF5L1, Snr1

Observed MW: 44 kDa Uniprot ID: Q12824

Application

Reactivity: Human

Tested Application: ELISA, IHC

Recommended dilution: IHC: 1:50-1:500

1


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>



Image:



Immunohistochemistry of paraffin-embedded human kidney tissue slide using FNab10186(SNF5 antibody) at dilution of 1:200