

anti- CST3 antibody

Product Information

Catalog No.:	FNab10197
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)
Rackground	

Background

None

Immunogen information

Immunogen:	cystatin C
Synonyms:	ARMD11, CST3, Cystatin 3, cystatin C, CYTC, Gamma trace, Post gamma globulin
Observed MW:	16 kDa
Uniprot ID :	P01034

Application

Reactivity:	Human, Mouse, Rat		
Tested Application:	ELISA, WB, IHC		
Recommended dilution: WB: 1:200-1:2000; IHC: 1:50-1:200			
Image:			


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)


Tel :(0086)027-87384275 Fax: (0086)027-87800889

www.fn-test.com





Immunohistochemistry of paraffin-embedded Rat testis tissue slide using FNab10197(CST3 antibody) at dilution of 1:200



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab10197(CST3 antibody) at dilution of 1:1000

2

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)