

# anti- SFRP2 antibody

### **Product Information**

Catalog No.: FNab10200

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

### **Background**

None

# Immunogen information

Immunogen: Secreted frizzled-related protein 2

Synonyms: FRP 2, FRP2, SARP 1, SARP1, SDF 5, sFRP 2, SFRP2,

UNQ361/PRO697

Observed MW: 33 kDa Uniprot ID: Q96HF1

## **Application**

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, WB, IHC

Recommended dilution: WB: 1:200-1:2000; IHC: 1:50-1:200

Image:


### 1

#### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using FNab10200(SFRP2 antibody) at dilution of 1:200



rat pancreas were subjected to SDS PAGE followed by western blot with FNab10200(SFRP2 antibody) at dilution of 1:1000