

# anti-SIR6 antibody

#### **Product Information**

Catalog No.: FNab10207

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

## **Background**

None

## **Immunogen information**

Immunogen: NAD-dependent protein deacetylase sirtuin-6

Synonyms: SIR2L6
Observed MW: 39 kDa
Uniprot ID: Q8N6T7

## **Application**

Reactivity: Human

Tested Application: ELISA, WB, IHC

Recommended dilution: WB: 1:200-1:2000; IHC: 1:50-1:200

Image:


#### 1

#### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Immunohistochemistry of paraffin-embedded human kidney cancer tissue slide using FNab10207(SIR6 antibody) at dilution of 1:200



Hela cells were subjected to SDS PAGE followed by western blot with FNab10207(SIR6 antibody) at dilution of 1:1000