

# anti- IBP2 antibody

#### **Product Information**

Catalog No.:	FNab10211
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	$\geq$ 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)
Background	

None

## **Immunogen information**

Immunogen:	Insulin-like growth factor-binding protein 2
Synonyms:	BP2, IBP 2, IBP2, IGF binding protein 2, IGF BP53, IGFBP 2, IGFBP2
Observed MW:	34 kDa
Uniprot ID :	P18065

### Application

Reactivity:	Human		
Tested Application:	ELISA, WB, IHC		
Recommended dilution: WB: 1:200-1:2000; IHC: 1:50-1:200			
Image:			


#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)


Tel :( 0086)027-87384275 Fax: (0086)027-87800889

www.fn-test.com





Immunohistochemistry of paraffin-embedded human kidney cancer tissue slide using FNab10211(IBP2 antibody) at dilution of 1:200



Hela cells were subjected to SDS PAGE followed by western blot with FNab10211(IBP2 antibody) at dilution of 1:1000

#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

www.fn-test.com