

anti- INLR1 antibody

Product Information

Catalog No.:	FNab10213
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	Interferon lambda receptor 1
Synonyms:	IL28RA, LICR2
Observed MW:	57 kDa
Uniprot ID :	Q8IU57

Application

Reactivity:	Mouse
Tested Application:	ELISA, IHC
Recommended dilution: IHC:	1:50-1:500
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded mouse heart tissue slide using FNab10213(INLR1 antibody) at dilution of 1:200