

anti- mapk11 antibody

Product Information

Catalog No.: FNab10226

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen: Mitogen-activated protein kinase 11

Synonyms: MAP kinase 11, MAP kinase p38 beta, MAPK 11, MAPK11, p38 2,

P38B, p38Beta, P38BETA2, PRKM11, SAPK2, SAPK2B

Observed MW: 41 kDa Uniprot ID: Q15759

Application

Reactivity: Human, Mouse

Tested Application: ELISA, WB, IHC, IF

Recommended dilution: WB: 1:500-1:2000; IHC: 1:50-1:500

Image:


1

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Immunohistochemistry of paraffin-embedded mouse brain tissue slide using FNab10226(mapk11 antibody) at dilution of 1:200




Immunofluorescent analysis of Hela cells using mapk11 antibody(FNab10226) at 1:50 dilution.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Jurkat cells were subjected to SDS PAGE followed by western blot with FNab10226(mapk11 antibody) at dilution of 1:1000