

# anti- akt1 antibody

### **Product Information**

Catalog No.: FNab10227

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

## **Background**

None

# Immunogen information

Immunogen: RAC-alpha serine/threonine-protein kinase

Synonyms: AKT, AKT1, AKT2, AKT3, PKB, PKB ALPHA, PRKBA, Protein

kinase B, Proto oncogene c Akt, RAC, RAC ALPHA, RAC PK alpha

Observed MW: 48/55 kDa
Uniprot ID: P31749

# **Application**

Reactivity: Human, rat

Tested Application: ELISA, WB, IHC

Recommended dilution: WB: 1:500-1:2000; IHC: 1:50-1:500

Image:


#### 1

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Immunohistochemistry of paraffin-embedded Rat kidney tissue slide using FNab10227(akt1 antibody) at dilution of 1:200



Hela cells were subjected to SDS PAGE followed by western blot with FNab10227(akt1 antibody) at dilution of 1:1000