

anti- mek2 antibody

Product Information

Catalog No.: FNab10241

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen: Dual specificity mitogen-activated protein kinase kinase 2

Synonyms: ERK activator kinase 2, MAP kinase kinase 2, MAP2K2, MAPK/ERK

kinase 2, MAPKK 2, MAPKK2, MEK 2, MEK2, MKK2, PRKMK2

Observed MW: 44 kDa Uniprot ID: P36507

Application

Reactivity: Human, Mouse, Rat Tested Application: ELISA, WB, IHC, IF

Recommended dilution: WB: 1:500-1:2000; IHC: 1:50-1:200; IF: 1:10-1:100

Image:


1

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Immunohistochemistry of paraffin-embedded rat kidney tissue slide using FNab10241(MEK2 antibody) at dilution of 1:100



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab10241(mek2 antibody) at dilution of 1:1000