

anti- HSP60 antibody

Product Information

Catalog No.: FNab10250

Size: $100\mu g$ Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

HSPD1, also known as HSP60, belongs to the chaperonin family and acts as a chaperone to enhance cell survival under physiological stresses. Hsp60 has been shown to be connected with many aspects of cell functions such as protein folding and assembling of polypeptide chains in mitochondria. Recently it has been reported that HSP60 is associated with apoptosis or inhibition of cancer cell growth.(21822415)

Immunogen information

Immunogen: heat shock 60kDa protein 1(chaperonin)

Synonyms: Heat shock protein 60, HSP 60, HSP60, HSP65, HSPD1

Observed MW: 60 kDa, 55 kDa

Uniprot ID: P10809

Application

Reactivity: Human, Mouse, Rat Tested Application: ELISA, WB, IHC

Recommended dilution: WB: 1:500 - 1:2000; IHC: 1:50 - 1:200

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>


1



Image:



Immunohistochemistry of paraffin-embedded rat kidney using FNab10250(HSP60 antibody) at dilution of 1:100



HeLa cells were subjected to SDS PAGE followed by western blot with FNab10250(HSP60 antibody) at dilution of 1:1000