

anti- OAT1 antibody

Product Information

Catalog No.: FNab10269

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24

months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen: solute carrier family 22 member 6

Synonyms: SLC22A6, HOAT1, OAT1, PAHT, ROAT1, solute carrier family 22

member 6

Observed MW: 60 kDa Uniprot ID: Q4U2R8

Application

Reactivity: Human, Mouse

Tested Application: ELISA, WB

Recommended dilution: WB: 1:200-1:2000

Image:


1

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Mouse kidney tissue were subjected to SDS PAGE followed by western blot with FNab10269(OAT1 antibody) at dilution of 1:1000