

anti- FYN antibody

Product Information

Catalog No.: FNab10273

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen: FYN oncogene related to SRC, FGR, YES

Synonyms: FYN, p59 Fyn, Proto oncogene c Fyn, Proto oncogene Syn, SLK, Src

like kinase, SYN, Tyrosine protein kinase Fyn

Observed MW: 55 kDa Uniprot ID: P06241

Application

Reactivity: Human

Tested Application: ELISA, WB, IHC

Recommended dilution: WB: 1:500 - 1:2000; IHC: 1:50 - 1:200

Image:


1

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab10273 (FYN Antibody) at dilution of 1:100



HeLa cells were subjected to SDS PAGE followed by western blot with FNab10273 (FYN antibody) at dilution of 1:1000