

anti- WDR74 antibody

Product Information

Catalog No.:	FNab10301	
Size:	100µg	
Form:	liquid	
Purification:	Immunogen affinity purified	
Purity:	\geq 95% as determined by SDS-PAGE	
Host:	Rabbit	
Clonality:	polyclonal	
IsoType:	IgG	
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 12 months (Avoid repeated freeze / thaw cycles.)	

Immunogen information

Immunogen:	WD repeat domain 74
Synonyms:	NOP seven associated protein 1, NSA1, WD repeat domain 74, WDR74
Observed MW:	40 kDa
Uniprot ID :	Q6RFH5

Application


Reactivity:	Human			
Tested Application:	ELISA, WB			
Recommended dilution:WB: 1:200 - 1:1000				
Image:				

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>