

# anti-PDCD10 antibody

## **Product Information**

Catalog No.: FNab10302

Size: 100μg Form: liquid

Purification: Protein A+G purification

Purity: ≥95% as determined by SDS-PAGE

Host: Mouse

Clonality: monoclonal

Clone ID: 2A8
IsoType: IgG1

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

## **Immunogen information**

Immunogen: programmed cell death 10

Synonyms: CCM3, PDCD10, programmed cell death 10, TFAR15

Observed MW: 25-30KD Uniprot ID: Q9BUL8

# **Application**

Reactivity: Human

Tested Application: ELISA, WB, IHC

Recommended dilution: WB: 1:500-1:2000; IHC: 1:50-1:200

Image:


#### 1

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Immunohistochemistry of paraffin-embedded human liver cancer using FNab10302(PDCD10 antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab10302(PDCD10 antibody) at dilution of 1:1000