

anti-TNFRSF6B antibody

Product Information

Catalog No.: FNab10343

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

Decoy receptor that can neutralize the cytotoxic ligands TNFS14/LIGHT, TNFSF15 and TNFSF6/FASL. Protects against apoptosis.

Immunogen information

Immunogen: tumor necrosis factor receptor superfamily, member 6b, decoy

Synonyms: DCR3, Decoy receptor 3, Decoy receptor for Fas ligand, M68,

TNFRSF6B, TR6, TNF6B

Observed MW: None
Uniprot ID: 095407

Application

Reactivity: Human

Tested Application: ELISA, IHC

Recommended dilution: IHC: 1:50-1:500

Image:


1

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Immunohistochemistry of paraffin-embedded human colon cancer using FNab10343(TNFRSF6B antibody) at dilution of 1:100

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>